


データ可視化ツール・IGVの紹介・実習



The screenshot shows the IGV website homepage. On the left is a navigation menu with links to Home, Downloads, Documents, Hosted Genomes, FAQ, IGV User Guide, File Formats, Release Notes, IGV for iPad, Credits, and Contact. Below the menu is a search bar and the Broad Institute logo. The main content area features a large banner with the IGV logo and a visualization of genomic data. Below the banner are sections for 'What's New' (dated September 2014), 'Overview' (describing IGV as a high-performance tool for exploring genomic data), 'Downloads' (with a registration link), 'Citing IGV' (providing citation information), and 'Funding' (listing funding sources like the National Cancer Institute). At the bottom, there are logos for various institutions and the URL <https://www.broadinstitute.org/igv/>.

なぜIGVを取り上げるか

データ可視化ツール

- ・自分のパソコン(ローカル環境)にインストールして使うタイプ
- ・サーバーに構築して、ネットワークで使うタイプ


The JBrowse Genome Browser

JBrowse is a fast, embeddable genome browser built completely with JavaScript and HTML5, with optional run-once data formatting tools written in Perl.

Featured Post

[Exploring structural variation using JBrowse](#) by Richard Finkers

Latest Release – JBrowse 1.11.6



コミュニティに広く利用、あるいはウェブ公開を目的とするには良いが、ネットワーク・情報セキュリティの高度な知識も要求される。

より大容量なデータに対応できる。

管理者的な人がいて、その人がやってくれるなら、これも良いが。

もっとお手軽なものとしてIGVを紹介

可視化ツールに求められるものは何か

膨大なデータを如何に直感的に理解できるようにするか
sortや絞り込みができる表データと対比双壁

- ・配列、GC ratio、遺伝子情報
- ・遺伝子発現情報
- ・SNPの位置情報・頻度情報
- ・様々なデータの精度情報

レファレンス配列 / gene model / gene annotationとNGSデータを並べて比較
複数のデータセットを並べて比較

色々なデータ(variant, 発現, ChIP, BSseq等々)を、様々なスケールで
比較・統合的に解釈できるようにしたい

ゲノムviewerに自分のデータを乗せ、
統合的直感的に解釈できること

可視化ツールをどう選ぶか

選択の基準

genome data viewing に求められるもの
取捨選択の基準

1. 無料 / 有料 / 基本無料
2. 個人的レベルの使用 / コミュニティーレベルの使用
3. 見るだけ/自分から色々工夫
4. アクセスのしやすさ・使いやすさ
導入に必要なコンピュータスペック
マニュアルは分かりやすいか
情報の多さ
利用の簡便さ
使っている人が近くにいるか

Integrative Genomics Viewer(IGV)

お手軽ツール

- ・アカデミックウェアで無料
- ・コミュニティーでの利用者が多いから、情報も多い
- ・javaのプログラムなので、オールプラットフォーム対応
- ・マニュアルは親切、サンプルデータのある
- ・WEBサーバーではなく、PCレベルでできる
- ・データ閲覧環境の共有が可能

誰もが簡便に使えるものが良い。

The screenshot shows the IGV website homepage. The sidebar on the left contains navigation links: Home, Downloads, Documents, Hosted Genomes, FAQ, IGV User Guide, File Formats, Release Notes, IGV for iPad, Credits, and Contact. Below the sidebar is a search bar and the Broad Institute logo. The main content area has a large banner with the text 'Integrative Genomics Viewer' and an image of the IGV interface. Below the banner are sections for 'What's New', 'Overview', 'Downloads', 'Citing IGV', and 'Funding'. The 'Citing IGV' section is highlighted with a red box and contains the following text:

Citing IGV

To cite your use of IGV in your publication:

James T. Robinson, Helga Thorvaldsdóttir, Wendy Winckler, Mitchell Guttman, Eric S. Lander, Gad Getz, Jill P. Mesirov. [Integrative Genomics Viewer \(IGV\): high-performance genomics data visualization and exploration](#). *Briefings in Bioinformatics* 14, 178-192 (2013).

The 'Funding' section mentions that the development of IGV is made possible by funding from the National Cancer Institute, the National Institute of General Medical Sciences of the National Institutes of Health, and the Starr Cancer Consortium. It also mentions that IGV participates in the GenomeSpace initiative, which is funded by the National Human Genome Research Institute. Logos for the National Cancer Institute, National Human Genome Research Institute, and GenomeSpace are shown at the bottom.

Integrative genomics viewer

James T Robinson, Helga Thorvaldsdóttir, Wendy Winckler, Mitchell Guttman, Eric S Lander, Gad Getz & Jill P Mesirov

[Affiliations](#) | [Corresponding authors](#)

Nature Biotechnology **29**, 24–26 (2011) | doi:10.1038/nbt.1754
Published online 10 January 2011

To the Editor:

Rapid improvements in sequencing and array-based platforms are resulting in a flood of diverse genome-wide data, including data from exome and whole-genome sequencing, epigenetic surveys, expression profiling of coding and noncoding RNAs, single nucleotide polymorphism (SNP) and copy number profiling, and functional assays. Analysis of these large, diverse data sets holds the promise of a more comprehensive understanding of the genome and its relation to human disease. Experienced and knowledgeable human review is an essential component of this process, complementing computational approaches. This calls for efficient and intuitive visualization tools able to scale to very large data sets and to flexibly integrate multiple data types, including clinical data. However, the sheer volume and scope of data pose a significant challenge to the development of such tools.

[Journal home](#)
[Current issue](#)
[For authors](#)
[Subscribe](#)
[E-alert sign up](#)
[RSS feed](#)


Subscribe today, save **50%** and receive **51** weekly issues of *Nature* in **print, online** and **mobile**.

Citations to this article

[Crossref \(10\)](#) [Scopus \(12\)](#) [Web of Science \(0\)](#)

Science jobs from [naturejobs](#)

Faculty Position

Harvard Medical School

Ramalingaswami Re-Entry Fellowship

Ministry of Science & Technology, Government of India



- [Home](#)
- [Downloads](#)
- [Documents](#)
- [Hosted Genomes](#)
- [FAQ](#)
- [IGV User Guide](#)
- [File Formats](#)
- [Release Notes](#)
- [IGV for iPad](#)
- [Credits](#)
- [Contact](#)

Search website

[Broad Home](#)
[Cancer Program](#)



© 2013 Broad Institute

Home



What's New



September 2014. The IGV iPad app can now be installed from the Apple App Store. *IGV for iPad* is a lightweight genomic data viewer that provides some of the functionality available in our regular desktop IGV. See the [IGV for iPad documentation](#) for details.

Overview



The **Integrative Genomics Viewer (IGV)** is a high-performance visualization tool for interactive exploration of large, integrated genomic datasets. It supports a wide variety of data types, including array-based and next-generation sequence data, and genomic annotations.

Downloads



Please [register](#) to download IGV. After registering, you can log in at any time using your email address. Permission to use IGV is granted under the GNU [LGPL license](#).

Citing IGV

To cite your use of IGV in your publication:

James T. Robinson, Helga Thorvaldsdóttir, Wendy Winckler, Mitchell Guttman, Eric S. Lander, Gad Getz, Jill P. Mesirov. [Integrative Genomics Viewer \(IGV\): high-performance genomics data visualization and exploration](#). *Briefings in Bioinformatics* **14**, 178–192 (2013).

Helga Thorvaldsdóttir, James T. Robinson, Jill P. Mesirov. [Integrative Genomics Viewer \(IGV\): high-performance genomics data visualization and exploration](#). *Briefings in Bioinformatics* **14**, 178–192 (2013).

Funding

Development of IGV is made possible by funding from the [National Cancer Institute](#), the [National Institute of General Medical Sciences](#) of the [National Institutes of Health](#), and the [Starr Cancer Consortium](#).

IGV participates in the [GenomeSpace](#) initiative, which is funded by the [National Human Genome Research Institute](#).



IGV User Guide


This guide describes the Integrative Genomics Viewer (IGV).

- To start IGV, go to the IGV downloads page: <http://www.broadinstitute.org/igv/download>.

[Look at a printer-friendly HTML version of the whole User Guide.](#)

- [User Interface](#)
- [Starting IGV](#)
- [Navigating](#)
- [Loading a Genome](#)
- [Viewing the Reference Genome](#)
- [Loading Data and Attributes](#)
- [Viewing Data](#)
- [Viewing Alignments](#)
- [Viewing Variants](#)
- [Gene List View](#)
- [Regions of Interest](#)
- [Sample Attributes](#)
- [Sorting, Grouping, and Filtering](#)
- [Saving and Restoring Sessions](#)
- [Server Configuration](#)
- [External Control of IGV](#)
- [Motif Finder](#)
- [igvtools](#)

[User Interface »](#)



igv Integrative Genomics Viewer

- Home
- Downloads**
- Documents
 - Hosted Genomes
 - FAQ
 - IGV User Guide
 - File Formats
 - Release Notes
 - IGV for iPad
 - Credits
- Contact

Search website

search


[Broad Home](#)
[Cancer Program](#)

BROAD
INSTITUTE

© 2013 Broad Institute

Home

Integrative Genomics Viewer



What's New



September 2014. The IGV iPad app can now be installed from the Apple App Store. *IGV for iPad* is a lightweight genomic data viewer that provides some of the functionality available in our regular desktop IGV. See the [IGV for iPad documentation](#) for details.

Overview



The **Integrative Genomics Viewer (IGV)** is a high-performance visualization tool for interactive exploration of large, integrated genomic datasets. It supports a wide variety of data types, including array-based and next-generation sequence data, and genomic annotations.

Downloads



Please [register](#) to download IGV. After registering, you can log in at any time using your email address. Permission to use IGV is granted under the GNU [LGPL license](#).

Citing IGV

To cite your use of IGV in your publication:

James T. Robinson, Helga Thorvaldsdóttir, Wendy Winckler, Mitchell Guttman, Eric S. Lander, Gad Getz, Jill P. Mesirov. [Integrative Genomics Viewer](#). *Nature Biotechnology* 29, 24–26 (2011)

Helga Thorvaldsdóttir, James T. Robinson, Jill P. Mesirov. [Integrative Genomics Viewer \(IGV\): high-performance genomics data visualization and exploration](#). *Briefings in Bioinformatics* 14, 178–192 (2013).

Funding

Development of IGV is made possible by funding from the [National Cancer Institute](#), the [National Institute of General Medical Sciences](#) of the [National Institutes of Health](#), and the [Starr Cancer Consortium](#).

IGV participates in the [GenomeSpace](#) initiative, which is funded by the [National Human Genome Research Institute](#).



レジストレーションが必要

Downloads

Integrative Genomics Viewer (IGV) (Version 2.3)

Install IGV

Options for installing and running IGV:

- (Mac only) Download and run the Mac application; or
- (All systems) Use the Java Web Start buttons (Mac users: see below for limitations); or
- (All systems) Download the binary distribution and run IGV from the command line.

1. Mac Application

Download and unzip the Mac App archive, then double-click the IGV application to run it. The application can be moved to the "Applications" folder, or anywhere else. **Note: This requires Java 7. Mac users with Java 6 (JRE 1.6) should use the binary distribution archive or the Java Web Start buttons below.**

Download
Mac App

2. Java Web Start

The buttons below use Java Web Start (JWS) to install and launch IGV directly from our web site.

***Mac Users:** The Java Web Start option is not recommended for Mac OSX Mountain Lion or higher. Using it requires that you set Gatekeeper security to its lowest level, and it is possible that even this will not be enough.

Chrome: Chrome does not automatically launch the Java Webstart files by default. Instead, the launch buttons below will download a "jnlp" file. This should appear in the lower left corner of the browser. Double-click the downloaded file to run.

Windows users: To run with more than 1.2 GB of memory you must install 64-bit Java. **Most Windows installs do not include 64-bit Java by default, even if the operating system is 64-bit.** Attempting to use the 2GB or greater launch options with 32-bit Java will result in the error "could not create virtual machine".

| | | | |
|------------------------------|----------------------------------------------------------------------------------------|-----------------------------------------------------------------------|----------------------------------------------------------------------------|
| Launch Launch with 750 MB | Launch Launch with 1.2 GB Maximum usable memory for Windows OS with 32-bit Java. | Launch Launch with 2 GB Maximum usable memory for 32-bit MacOS. | Launch Launch with 10 GB For large memory machines with 64-bit Java. |
|------------------------------|----------------------------------------------------------------------------------------|-----------------------------------------------------------------------|----------------------------------------------------------------------------|

3. Binary Distribution


Download and unzip the binary distribution archive in a folder of your choosing. IGV is launched from a command prompt – follow instructions in the "readme" file. To launch igv on Mac or Linux platforms use the shell script "igv.sh". On Windows use "igv.bat".

Download
Binary Distribution


igvtools

Utilities for preprocessing data files.

- [igvtools_2.3.40.zip](#)



いくつかの生物種は登録されている



新たな生物種をロードすることも可能

ゲノムViewerなので次世代DNAシーケンサーのデータに限定されない。
マイクロアレイの結果や、ゲノムアノテーションの情報も随時表示できる。

対応するファイル形式に応じて、表示方法が決まる。

File Formats


- File Extension Identifies Format
- Recommended File Formats
- BAM
- BED
- BedGraph
- bigBed
- bigWig
- Birdsuite Files
- broadPeak
- CBS
- CN
- Custom File Formats
- Cytoband
- FASTA
- GCT
- genePred
- GFF/GTF
- GISTIC
- Goby
- GWAS
- IGV
- LOH
- MAF (Multiple Alignment Format)
- MAF (Mutation Annotation Format)
- Merged BAM File
- MUT
- narrowPeak
- PSL
- RES
- SAM
- Sample Information
- SEG
- SNP
- TAB
- TDF
- Track Line
- Type Line
- VCF
- WIG
- chrom.sizes

File Formats

IGV supports a number of different file formats for experimental data and genome annotations. For a complete list of supported formats see <http://www.broadinstitute.org/igv/FileFormats>. The following table shows the recommended file formats for a number of common data types.

| Source Data | Recommended File Formats |
|-------------------------|----------------------------------------|
| ChIP-Seq, RNA-Seq | WIG, TDF |
| Copy number | CN, SNP, TDF, canary_calls (Birdsuite) |
| Gene expression data | GCT, RES, TDF |
| Genome annotations | GFF, BED, GTF, PSL, UCSC table format |
| GISTIC data | GISTIC |
| LOH data | LOH, TDF |
| Mutation data | MUT, MAF |
| Variant calls | VCF |
| RNAi data | GCT |
| Segmented data | SEG, CBS |
| Sequence alignment data | BAM, SAM, PSL |
| Any numeric data | IGV, WIG, TDF |
| Sample metadata | Tab-delimited sample info file |

公開情報のviewerとして



その他の便利機能

セッションの保存

表示しているデータの読み込み状況を、それごと保存。

セッションをロードすることで、意図した画面を表示できる。

データセットが揃っていること、フォルダー構造が同一である必要がある。


バッチ処理

重要領域の画面スナップショットを自動で取ったりできる。


```
new
load myfile.bam
snapshotDirectory mySnapshotDirectory
genome hg18
goto chr1:65,289,335-65,309,335
sort position
collapse
snapshot
goto chr1:113,144,120-113,164,120
sort base
collapse
snapshot
```

Inversions


An inversion is a large section of DNA that is reversed in the subject genome compared to the reference genome.



When an inversion shows up in paired-end reads, the reads are distinctively variant from the reference genome.



This appears in IGV as shown below.






Interpreting Color by Insert Size

The inferred insert size can be used to detect structural variants, such as:


- deletions
- insertions
- inter-chromosomal rearrangements

IGV uses color coding to flag anomalous insert sizes. When you select Color alignments>by insert size in the popup menu, the default coloring scheme is:


-  for an insert that is larger than expected
-  for an insert that is smaller than expected
-  for paired end reads that are coded by the chromosome on which their mates can be found

Deletions


A deletion is a large section of DNA that is absent in the subject genome compared to the reference genome.



The "expected" insert size is the insert size obtained in sequencing the subject genome. The "inferred" insert size is the insert size that would result in the reference genome, assuming the same pair of reads.




RNA-Seqのデータ表示させる





GCTファイルでgene ローカスの発現情報を図示


| # | | | | |
|---|-----------|---------------------|-----------|-----------|
| # | Name | Description | 2D_2D2L | 4D_2D2L |
| | ANAC001 | @Chr1:3630-5899 | -2.60184 | -2.60956 |
| | DCL1 | @Chr1:23145-33153 | -0.742675 | -1.5642 |
| | MIR838A | @Chr1:23145-33153 | 0 | 0 |
| | AT1G01073 | @Chr1:44676-44787 | 0 | 0 |
| | IQD18 | @Chr1:52238-54692 | -1.93871 | -1.13128 |
| | AT1G01115 | @Chr1:56623-56740 | 0 | 0 |
| | GIF2 | @Chr1:72338-74737 | -0.251287 | -0.616679 |
| | AT1G01180 | @Chr1:75582-76758 | 0.45929 | -0.809567 |
| | AT1G01210 | @Chr1:88897-89745 | 1.6964 | 0.857196 |
| | FKGP | @Chr1:91375-95651 | -0.174589 | 0.725947 |
| | AT1G01240 | @Chr1:99893-101834 | -0.226384 | -0.936641 |
| | AT1G01260 | @Chr1:108945-111609 | -0.161848 | 0.315699 |
| | CYP703A2 | @Chr1:112262-113947 | 0 | 0 |
| | CNX3 | @Chr1:114285-116108 | 0.111249 | -0.551359 |
| | AT1G01300 | @Chr1:116942-118764 | -0.68348 | 0.108578 |



Gene listを定義して
サンプルごと
条件ごと
の発現・発現変動を
カラーマップできる




IGV実習



IGVの使用法を学ぶと共に
先のファイルフォーマットも
確認しよう

以下のファイルを確認

```
buc.genome.fasta  
buc.gtf  
buc_cg.wig  
illumina_ex_B2_Read_bowtie2.mate.sort.bam  
illumina_ex_B2_Read_bowtie2.mate.sort.bam.bai  
illumina_ex_B4_Read_bowtie2.mate.sort.bam  
illumina_ex_B4_Read_bowtie2.mate.sort.bam.bai
```



登録されていない生物種・配列でも、自分でimportすればOK


Buchneraゲノムを使う

Genome -> Load Genome from File

buc.genome.fasta


を読み込む

ビューイングの拡大縮小



カーソルのgene modelに
持っていくことで、gtfファイル
の記述がpop upされる。



トラックの移動
マウスのドラッグ移動が可能
← →による移動も可能
Feature trackを指定して、
Ctrl+F gene model単位で右に移動
Ctrl+B gene model単位で左に移動
Shift+Ctrl+F exon単位で右に移動



シーケンスデータを読み込ませる
File -> Load from File

illumina_ex_B2_Read_bowtie2.mate.sort.bam

同ファイル名のindexファイル
illumina_ex_B2_Read_bowtie2.mate.sort.bam.bai
が同じフォルダーにないといけない。



こういう所を
クリックしたまま
移動させることで
幅を変更

WIGファイルを読み込む
 ウィンドウサイズ50baseでCG ratioを計算したWIGファイルを作成してある
 buc_cg.wig

```

variableStep chrom=gn:buc
25      16
26      16
27      16
28      16
29      16
30      14
31      12
32      12
33      10
34      10
35      8
36      8
37      8
38      8
  
```

File -> Load from File でそのまま読み込ませることも可能だが、
 テキストファイルで膨大なサイズだと非常に時間がかかる。
 Tools -> Run igvtools

IGVtoolsが起動

Command Title

Input File G:\training_course\2013春\配布_Data\05_NGS_IGV\buc_cg.wig

Output File G:\training_course\2013春\配布_Data\05_NGS_IGV\buc_cg.wig.tdf

Genome G:\training_course\2013春\配布_Data\05_NGS_IGV\buc_genome.fasta

Tile and Count Options

Zoom Levels ?

Window Functions ☐ Min ☐ Max ☒ Mean ☐ Median

☐ 2% ☐ 10% ☐ 90% ☐ 98%

Probe to Loci Mapping

Window Size 25



Sort Options



Temp Directory



Max Records 500000


Close Run

Messages









IGV紹介のまとめ

可視化ツールとして十分な機能を持つ

- ・無料
- ・比較的簡単・お手軽
- ・自分で見るとにもよし、人に見せるためにもよし
- ・利用範囲は次世代DNAシーケンサーに限定しない
広くゲノミクスの解析に有用

ごく一部のみの機能を紹介しました。
ウェブサイトを見ながら復習をお勧めします。